



Software Feature	Medilog Darwin full version	Medilog Darwin ScanLab	Medilog Darwin Lite (Network Licence)
ADAPT™	99.9% accuracy	99.9% accuracy	99,9% accuracy
Workflows	Drop down menu for workflow selection	Workflow is fixed. Screen layouts can be changed.	Depending on which Darwin version it is connected to
Upload & Import Possibilities	Yes	Yes	No
ST Analysis	Yes	Yes	Yes
QT Analysis	Yes	Yes	Yes
HRV Analysis	Time Domain Frequency Domain	Time Domain	Depending which Darwin version it is connected to
HRV-Spectrogram “Fire of Life™”	Yes	No	
Pacemaker Analysis	Yes	Yes	Yes
Atrial Analysis	Yes	No	Depending which Darwin version it is connected to
Apnoea Analysis	Yes	No	
ECHO View® Module	Yes	No	
Arrhythmia Overview Module	Yes	No	
Time Setup	Yes	No	
Reports	Full Set	No reports for Atrial Analysis, Apnoea Analysis, HRV Frequency Domain	
HL7 interface	Option	Option	Option
SEMA Connectivity	Yes	Yes	Yes
Export possibility via ISHNE, MatLab or MIT	Yes	Yes	Yes