



SCHILLER
The Art of Diagnostics

Press Information

SCHILLER CARDIOVIT CS-200 – new design and additional features

Baer, July 2009 Passion fuels innovation! A passion for pioneering inventions has always distinguished the products of SCHILLER's cardiopulmonary diagnostic line. SCHILLER proudly introduces its redesigned complete diagnostic solution, the CS-200.

7 functions available in one system:

- Resting ECG including measurement and serial test comparison
- Exercise testing with ergometers and treadmills
- 24-hour ambulatory ECG (Holter)
- 24-hour ambulatory BP (Holter)
- Spirometry PFT
- On-board data management SEMA-200
- Networking capabilities



Discreet but important external modifications encompass a large, swivel-mounted 19" monitor, as well as an ergonomic design.

A new, simplified user-interface as well as various new holders and mounting kits for external devices such as gel bottles, bar code scanners, spirometry sensors and defibrillators provide even more comfort for both doctor and patient.

**For more information, please visit our website:
<http://www.schiller.ch>**